


Material / Materiaal			
No./Nr.	Plan	Article / Artikel	Number / Aantal
1	Y001468	1200701	1
		1200754	1
2	Y001440	1200601	1
		1200654	1
		1001448	1
3	---	1001448	1
4	X001016	1001682	1
5	Y000058	F000120	2
6	---	2060000	2
6A	X000578	1000967	2
6B	Y000202	1059900	2
6C	Y000199	1059800	2
6D	∅3x16 DIN1481	1003791	2
6E	∅3x20	1999776	2
6F	X001240	1004551	2
7	M5x30 DIN7985	1001475	4
8	M5 DIN125	1000972	8
9	M5 DIN985	1001679	4
10	X001069	1001633	2

WHEEL CHAIR 101 CROSS 50#				Plan Number
				EC00487/A
Scale: 1:1	Drafted	Verified	Confirmed	
Date	03/06/2020	03/06/2020	03/06/2020	
Initials	AL	JJQ	TG	