



Material			
No.	Plan	Article	Quantity
1	(YP06091 ILLICO)	2794182	1
	(YP06091 VITA BED)	2794113	1
	(YP06091 LUNA2)	2794137	1
2	XP01909	2993605	1
3	XP01913	2993620	1
4	Hand Grip	1906253	1
5	Buckle	1906981	1
6	Belt	2993608	1
7	Belt	2993609	1
8	Tube plug 42,4	1911703	1
9	X003559	1998689	2
10*	Sticker	T001035	1
11**	YP05298	2706601	1

* - not depicted / nie ujeta

** - option / opcja

EP00268 / F

SELF RAISING AID FOR ILLICO / VITA BED / LUNA 2				Plan No.
				EP00268/ F
Scale: 1:6	Drafted	Verified	Confirmed	
Date	09/11/2017	10/09/2018	10/09/2018	
Initials	PM	AZ	TK	