

Detail A
Scale: 1:3


EP00198 / G

FOOT SIDE 997 (Luna 2 Luna Basic 2)				Plan No.
				1/2 EP00198/ G
Scale: 1:6	Drafted	Verified	Confirmed	
Date	26/03/2015	02/03/2016	02/03/2016	
Initials	AS	MKR	TK	



Material			
No.	Plan	Article	Quantity
1	YP05064 DP	2899037	1
2	YP05106	2899637	1
3	YP05078	2899737	1
4	---	1909002	2
5	Linak motor 50	1996202	1
6	---	1998885	1
7	M8x16	2990074	2
8	Cup20x20	1001612	4
9	M8X25 DIN 933	1001442	2
10	X002480	1993793	4
11	X001016	1001682	8
12	M8x16 DIN933	1900449	2
13	8.4X17X1.6	1003199	4
14	M8x20 DIN933	1009400	2
15	4.8X19 DIN7983	1001534	2
16	Y012035	H810065	1
17	D10X40	1991481	1

EP00198 / G

FOOT SIDE 997 (Luna 2 Luna Basic 2)				Plan No.
2/2				EP00198/ G
Scale: 1:6	Drafted	Verified	Confirmed	
Date	26/03/2015	02/03/2016	02/03/2016	
Initials	AS	MKR	TK	