



Material			
No.	Plan	Article	Quantity
1	X002699	1996164 (L)	1
	X002853	1996165 (R)	1
2	M6x40 ISO7380	1006111	4
3	Y009175	1641654	3
4	M6 DIN125	1001564	8
5	M6 DIN125	1000973	8
6	M6 DIN985	1001680	4
7	M8x16 ISO7380	1991858	4
8	Y009176	1641754 (R)	1
	Y011973	1660154 (L)	1
9	M6x10	1003658	1
10	M10 DIN934	1001672	1
11	M10x19	1991432	1
12	M6x50 ISO7380	1006997	4
13	YP01682	1876454 (R)	1
	YP01683	1876554 (L)	1
14	XP00411	2990439	2
15	---	1009905	4
16	---	1001612	4
17	Y003626	1476554	1
18	Y013589	1681054 (R)	1
	Y013590	1681154 (L)	1
19	M5x16 ISO7380	1006112	4
20	Y003955	1496000	2
21	Y004459	1513354 (R)	1
	Y004458	1513254 (L)	1
22	M8x15	1006661	2
23	---	1903921	1
24*	---	1909020	1
25*	Y019535	1785854 (R)	1
	Y019593	1785954 (L)	1
26*	M5x10 ISO7380	1007063	2

* For LE drive only

E000480 / E

ARMREST RAPIDO / EXPRESS 2000				Plan No.
				E000480 / E
Scale : 1:4	Drafted	Verified	Confirmed	
Date	18/07/2008	07/02/2013	07/02/2013	
Init.	AB	SD	PJ	