



Detail B
Scale: 1:3

EP00211 / D

HEAD SIDE WITH WOODEN LATHS				Plan No.
LUNA 2 (997)				EP00211 / D
Scale: 1:6	Drafted	Verified	Confirmed	1/2
Date	08/01/2016	14/09/2020	14/09/2020	
Initials	KR	MKR	TK	



Material			
No.	Plan	Article	Quantity
1	YP05077_DP	2899237	1
2	YP05229	2704837	1
3	YP05110	2899800	2
4	X001016	1001682	2
5	M8x45_DIN933	1005901	2
6	XP01421	2993116	7
7	XP01420	2993117	14
8	---	1001637	2
9	D8,4 DIN125	1003199	4
10	---	1998885	1
11	---	1991481	1
12	M8X16	1991858	4
13	---	1902971	1
14	---	1995837	1
15	---	1905445	1
16	Motor LA27	1996201	1
17	---	1996203	1
18	---	1996204	1
19	---	1996205	2
20	---	1005899	2
21	---	1996206	1
22	X003329	1996776	2
23*	---	W620228	1

* - not depicted / nie ujeta

EP00211 / D

HEAD SIDE WITH WOODEN LATHS			Plan No.
LUNA 2 (997)			EP00211 / D
Scale: 1:6	Drafted	Verified	Confirmed
Date	08/01/2016	14/09/2020	14/09/2020
Initials	KR	MKR	TK
			