


Bill of Material			
For 1) PH00 2) PH01 3) PH03 4) PH04			
No.	Plan	Article	Quantity
1	M6x35 ISO7380 (Black)	2994131	4/2*
2	M6 DIN985 (Black)	2994063	4
3	D6,4 DIN125 (Black)	2994414	4
PH00			
4	Cap D22	1990648	2
5	P00048919 H30	2669065	2
	P00049631 H35	2669165	2
	P00049619 H40	2669265	2
PH01			
6	---	2993810	2
7	D4,8x16 DIN7983	1001533	2
8	X001473	1002994	2
9	P00048818 H30	2669465	2
	P00049957 H35	2669565	2
	P00049627 H40	2669665	2
PH03			
7	D4,8x16 DIN7983	1001533	2
8	X001473	1002994	2
10	EP00246	2993853	2
11	M6x10 ISO7380 (Black)	2994127	4
12	P00049712 H30	2669365	2
	P00048919 H35	2669065	2
	P00049631 H40	2669165	2
PH04 (Backerst D)			
13	GRIP BLACK	1908236 (39-46)	2
	GRIP BLACK	3994411 (48-50)	2
14	P00023273	2638354	2
15	P00027229	2643800	2
16	D3x16 DIN1481	1003791	2
17	Sterknob M6x15	1906622	2
18	X004876	1747754	2
19	Y017422	1747954	2
PH05			
20	EP00354	2994210	1



BACKREST / V500 ACTIVE			
Scale: 1:6	Drafted	Verified	Confirmed
Date	2023-01-12	17/02/2023	17/02/2023
Initials	AS	AZ	AS

* - Replace by 2994086 (x2) for B82S

Plan No.
E00049624/B1.1