


No.	Plan	Article	Qty
1	P00033432	2995478	1
2	P00056219	2674837	4
3	---	---	---
4	M6x20 ISO7380	1006999	2
5	M6x12 ISO7380	1007000	10
6	P00058649	US00016	1
7*	P00058910	U000452	1
8*	P00058912	US00017	4
9*	P00062587	2997908	2



* - only for wooden panels ordered separately (without bed)

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP

CLUB VARIO D	Remove 2x1008492 and 2x2990295 before installing 2995478
CLUB VARIO	Refer to drawing E00045757

ASSEMBLY FRONT/REAR PANEL / CLUB VARIO				---
				Plan No. E00045758/D1.1
Scale: 1:8	Drafted	Verified	Confirmed	
Date	2022-10-28	05/02/2024	05/02/2024	
Initials	OY	AZ	AZ	