


Detail A

2991419

Sheet 1/2

LOWER FRAME INOVYS II 16 INCH				Plan No. E000713/I
ONDERFRAME INOVYS II 16 INCH				
Scale: 4:25	Drafted	Verified	Confirmed	
Date	27/07/2011	19/01/2023	19/01/2023	
Initials	SD	DR	AS	

E000713 / I

Material			
No.	Plan	Article	Quantity
1	X000188	G078500	2
2	X001016	1001682	1
3	X003366	1612600	2
4	X003407	1997399	2
5	X003451	1658800	2
6	X003452	1658900	2
7	X003596	1998528	2
8	Y004353	F000411	2
9	Y010736	16386.. (40)	1
	Y011927	16596.. (45)	1
	Y011934	16597.. (50)	1
	YP06450	27342.. (35)	1
10	Y011640	1655100	1
11	Y011652	1655200	1
12	Y017901	1758901	2
13	Y016422	1000181	2
14	M12 DIN125	1001561	2
15	---	1001610	2
16	M6 DIN985	1001680	3
17	M12 DIN985	1001684	4
18	M16x35 ISO7380	1006057	1
19	M6x40 ISO7380	1006111	2
20	M12x75 DIN931	1009991	2
21	M7	1900008	1
22	200x50	1900184	2
23	6001-2Z	1992559	4
24	---	1992599	2
25	M8x45 ISO7380	1994473	1
26	M8x76 DIN912	1998993	2
27	M8	1998994	2
28	400N	1005602	1
29	Left	1900099	1
30	Right	1900100	1
31	Left	1910125	1
32	Right	1910126	1
33	M12 DIN433	2993281	2

Material			
No.	Plan	Article	Quantity
34*	X004032	1901832	1
35*	X002565	1993963	1


34* - use for frame width 35 with seat height 42

35* - use for frame width 35 with seat height 39

Sheet 2/2

2991419

E000713 / I

LOWER FRAME INOVYS II 16 INCH ONDERFRAME INOVYS II 16 INCH				Plan No.
				E000713/I
Scale: 4:25	Drafted	Verified	Confirmed	
Date	27/07/2011	19/01/2023	19/01/2023	
Initials	SD	DR	AS	