



Material			
No.	Plan	Article	Quantity
1	X003872	1901134	1
2	Y017946	1759204	1
3	M6 DIN125A	1000973	2
4	M6x20	1006999	2
5	M6x12	1008492 *	2
	M6 DIN985	1001680 **	2
	M5 DIN985	1001679 ***	2
6	M6x60	1992915 *	2
	M6x10	1900641 **	2
	M5x16	1006112 ***	2

* for B02/B03 armrest
 ** for B04/B05 armrest
 *** for BI5 V500/B04 V300DL

- (1) Left B02/B03
- (2) Right B02/B03
- (3) Left B04/B05
- (4) Right B04/B05
- (5) Left BI5 V500/B04 V300DL
- (6) Right BI5 V500/B04 V300DL

***	**	**	**	*	*
1914428 (6)	1914427 (5)	1907411 (4)	1907410 (3)	1907409 (2)	1907408 (1)

E000739 / E

B66 TYPE V-SERIES ADJUSTABLE				Plan No.
				E000739 / E
Scale: 2:5	Drafted	Verified	Confirmed	
Date	19/07/2011	26/10/2017	26/10/2017	
Initials	TT	TT	PJ	