

Sheet 1/2

2990906

E000629 / J

<h3>ADJUSTABLE FOOTREST B06 V200/V300</h3>				Plan No. / Plannr.
				E000629 / J
Scale / Schaal: 3:10	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	28/04/2010	02/07/2020	02/07/2020	
Init. / Init.	TM	SD	PJ	

Material / Materiaal			
No. / Nr.	Plan	Article / Artikel	Quantity / Aantal
1	X004414	1904130	1
2	X004416	1904202	1
3	X004417	1719200	1
4	X004418	1904207	1
5	X004423	1717300	1
6	X004570	1904646 (L)	(39-42) 1
	X004571	1904647 (R)	(39-42) 1
	X004457	1904475 (L)	(44-46) 1
	X004453	1904210 (R)	(44-46) 1
	X004572	1904648 (L)	(48-60) 1
	X004573	1904649 (R)	(48-60) 1
6*	XC00158	3991003 (L)	(39) 1
	XC00157	3991002 (R)	(39) 1
7	X004464	1904212	1


* Only for Express

Material / Materiaal			
No. / Nr.	Plan	Article / Artikel	Quantity / Aantal
8	X004465	1904213	1
9	X004658	1904968	1
10	X004660	1905195	1
11	X004687	1905639	1
12	Y015007	17300..	1
13	Y016499	F006269	1
14	Y016501	1733500	1
15	---	1001680	2
16	M6x50 DIN912	1007970	1
17	M6x40 DIN912	1003909	2
18	---	1001494	1
19	---	1008492	1
20	---	1009588	1
21	---	1905640	1
22	---	1905642	1
23	---	1001564	1

Sheet 2/2

2990906

E000629 / J

ADJUSTABLE FOOTREST B06 V200/V300				Plan No. / Plannr.
				E000629 / J
Scale / Schaal: 3:10	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	28/04/2010	02/07/2020	02/07/2020	
Init. / Init.	TM	SD	PJ	