



EXPLODED VIEW 276				Plan No. / Plannr.
1/2				E000994C
Scale / Schaal: 1:8	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	09/07/2015	15/05/2017	15/05/2017	
Init. / Init.	NV	DB	TK	

Material			
No. / Nr.	Plan / Rysunek	Article / Artikel	Number / Aantal
1	Y004065	15071..	1
2	Y004075	15075..	1
3	Y004067	15072..	1
4	Y004072	15073..	2
5	Y004433	15076..	2
6	---	1993412	1
7A, 7B	Complete / Komlet	1995385 ---	2 ---
8	M12x60	1993419	2
9A, 9B, 9C	Complete / Komlet	1993416 ---	2 ---
10A, 10B	Complete / Komlet	1995386 ---	2 ---
11	X002308	1507700	2
12	X001016	1001682	2
13	M6x75 DIN 912/8.8	1993413	2
14	X002306	15023..	2
15	X002443	15171..	2
16	X002448	1517200	2
17	D6,4 DIN 125 A	1001564	6
18	---	1003558	2
19	M6 DIN985	1001680	10
20	M6x40	1006111	2
21	M6x60 DIN 912/8.8	1005957	4
22	---	1000973	4
23	X000106	G003218	4
24	D10,5 DIN 125 A	1001558	2
25	M6x16	1008459	2
26	X002305	1001545	4
27	X002305	1993120	2
28	---	H016007	2
29	---	1992600	2
30	X000586	1005742	1
31	---	1993417	1

Material			
No. / Nr.	Plan / Rysunek	Article / Artykul	Number / Iloscl
32	---	1992400	2
33	M12 DIN 127 A	1001572	2
34	---	1909113	1
35	XP00071	2990040	2
36	---	1990816	4
37	---	1003901	2
38	---	1994587	1
39	X002335	1993411	1
40	M5x10 DIN 916/45H	1001485	2
41	L=955 mm	1003279	2
	L=1020 mm	1004094	2
42	---	1996192	2
* 43	---	T001013	2
* 44	---	T001014	1
* 45	---	PP00381	1
* 46	LOGO REHA-POL-A	T001002	1
* 47	LOGO VERMEIREN	T001003	2
* 48	---	W620150	2
* 49	---	W620151	2
* 50	---	1003558	2

*Not pictured.

EXPLODED VIEW 276				Plan No. / Plannr.
				2/2 E000994 C
Scale / Schaal:	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	09/07/2015	15/05/2017	15/05/2017	
Init. / Init.	NV	DB	TK	